

ID	First Name	Last Name	Sex	Age	City	Phone
12546	George	Best	M	31	Chicago	254876
54275	Miriam	Mathew	F	25	Boston	378177
24868	Anandan	Namaka	M	15	Raleigh	146683
86736	Barumah	Joy	F	14	Chicago	257986
66629	Amanda	Bell	F	28	Chicago	257984
54546	Ana	Ivanka	F	26	Chicago	255455
86773	Joy	Mukherji	M	14	Chicago	257987
98768	Mala	Sinha	F	45	Chicago	258786
74668	Jasmine	Rose	F	37	Chicago	257987
45355	Mohanam	Manmohan	M	60	Chicago	257923
35465	Justin	Lake	M	10	Chicago	257487
78656	Manny	Fred	M	85	Chicago	254987
13549	Elton	George	M	21	Chicago	257836
87866	Veronica	Cooper	F	40	Chicago	257967
65466	Madge	Hoper	F	15	Chicago	253584
69978	Marekel	Noikle	M	45	Chicago	253248
43527	Invan	Vankek	M	55	Chicago	257646
24868	Anandan	Namaka	M	15	Chicago	146683
86736	Barumah	Joy	F	14	Chicago	257986
66629	Amanda	Bell	F	28	Chicago	257984
54546	Ana	Ivanka	F	26	Chicago	255455
86773	Joy	Mukherji	M	14	Chicago	257987
98768	Mala	Sinha	F	45	Chicago	258786
74668	Jasmine	Rose	F	37	Chicago	257987
45355	Mohanam	Manmohan	M	60	Chicago	257923
35465	Justin	Lake	M	10	Chicago	257487
78656	Manny	Fred	M	85	Chicago	254987
13549	Elton	George	M	21	Chicago	257836
87866	Veronica	Cooper	F	40	Chicago	257967
65466	Madge	Hoper	F	15	Chicago	253584
69978	Marekel	Noikle	M	45	Chicago	253248
43527	Invan	Vankek	M	55	Chicago	257646

## PATIENT INFORMATION MANAGER:

- Patient details includes patient ID, name, address, sex, age, contact info, etc
- Option to auto create Patient IDs
- Patient ID is customizable to each imaging facility
- Gender based validation of study data. Eg: User will not be allowed to select Obstetrics study for a male patient
- Sort patient info based patient ID or patient name
- Search patient based on ID, name, study type, date, referring physician, etc.

## STUDY MANAGER:




- Study details include Study UID, Study ID, Accession number, modality, referring physician, location, etc.
- Add multiple modalities to be connected to the worklist generator.
- Pre stored modality and study list.
- Add new study type based on modality name, sex and study name
- Allows only authorized users to access the information
- Schedule patient studies
- Location indicator for multiple modalities. Eg: CT in room1, CT in room 2, etc.
- Track patient visits and scans for different modalities

## ADMINISTRATION MODULE:

- Security logs for tracking errors
- Manage multiple DICOM devices

## COMPATIBILITY:

- Compatible with any DICOM modality that supports worklist query
- Compatible with DICOM converters for non-DICOM modalities, that support querying from a worklist server



## SYSTEM REQUIREMENTS:

- Operating System: Windows 2000 (SP4) or Windows XP (SP2)  
Processing Power: Pentium IV 2.0 Ghz or more.  
Memory: 512 MB RAM or more.  
Hard Disk Space: 50 MB of Free Disk Space for installation.



### CORPORATE OFFICE

ASHVA TECHNOLOGIES PVT. LTD.

#15, ARK Colony,  
'Urmilla House', - 1st Floor, Eldams Road, Alwarpet, Chennai -  
600 018. India.  
Tel : 91 - 44 - 4208 8867 / 76, Fax : 91- 44 - 42088897

### ASHVA TECHNOLOGIES LTD.

505 Hwy 169 N,  
Waterford Park, Suite 300,  
Plymouth, MN 55441, U.S.A  
Tel: 763 746 9223  
Fax: 763 746 9224  
Email: marketing@ashvatech.com

### ASHVA TECHNOLOGIES.

Sdn. Bhd. (742382-M)  
No. 104 A, Jalan SS 24/2,  
Taman Megah, 47301 Petaling Jaya,  
Selangor Darul Ehsan, Malaysia.  
H/P: (6017) 2340805

[www.ashva.com](http://www.ashva.com)